

I N S T R U C T I O N S

Flowmi™ Tip Strainers

Catalog No. H13680-0040 – 40µm Porosity
H13680-0070 – 70µm Porosity

Figure 1

ENGLISH

The Flowmi™ Tip Strainers are pre-sterilized, disposable strainers that fit most P1000* pipette tips and facilitate removal of large particles which may cause clogs during flow cytometry analyses.

Tip Strainer Tray and Packaging:

Flowmi™ Tip Strainers are individually arranged in a tray to enable direct press and fit attachment to your pipette tip. Open the package by tearing at the notched location. Slight downward pressure on the cover will loosen the catch and allow the cover to slide from side to side. *Note: Tip Strainers are sterile and if opened under laminar flow hood can be resealed in the zip-closure bag to maintain sterility.*

Directions for Use (see Figure 1):

- Using a P1000 pipettor, select the desired volume, affix a 1000µl pipette tip, and aspirate the sample.
- While the sample is still inside the pipette tip, push it into one of the Flowmi™ Tip Strainers and press down to secure the joining of the tip to the strainer. Use the same force to affix the Flowmi™ tip as you use to attach pipette tip to pipettor. You may want to practice one or two times with deionized water to be sure you know the correct force to attach Flowmi™
- Dispense the liquid from the pipette so it flows through the tip strainer and into your receptacle.** The sample has now been successfully filtered and is ready for analysis on a flow cytometer (or FACS provided these steps have all been performed in a laminar flow hood).
- Eject** the used pipette tip with the attached tip strainer into an appropriate waste receptacle.

Technical Specifications:

Flowmi™ Tip Strainers are made in a unique process that transforms solid polyethylene (PE) into a porous material with a precise pore size and a specific shape to fit most 1000µl pipette tips.*

- Pre-sterilized by gamma-irradiation for one-time use and disposal (Do Not Autoclave)
- Tip Strainers and tray are made from FDA grade materials
- Produced in the USA

*Flowmi™ Tip Strainers have been tested and fit the following pipette tips:

- Fisherbrand Sure-one
- VWR brand (53510-082)
- Gilson
- Axygen (TF-1000)
- Eppendorf epTIPS

If using a pipette tip other than shown above, we recommend you test fit prior to filtering a sample.

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Flowmi™ is a trademark of Bel-Art® – SP Scienceware®.
© 2018. All Rights Reserved.
Bel-Art® – SP Scienceware® is a division of SP Industries, Inc.