

INSTRUCTIONS



Switch-Grid Racks



CATALOG No.	TUBE SIZE	HALF RACK CAPACITY	FULL RACK CAPACITY	GRID COLOR
18745-0010	10-13mm	42 (7 x 6)	84 (14 x 6)	Orange
18745-0011	13-16mm	30 (6 x 5)	60 (12 x 5)	Blue
18745-0012	16-20mm	20 (5 x 4)	40 (10 x 4)	Green
18745-0013	20-25mm	12 (4 x 3)	24 (8 x 3)	Fuchsia
18745-0014	25-30mm	6 (3 x 2)	12 (6 x 2)	Yellow



Rack components:


- Switch-Grid Rack Frame
- Grids (4/pack)
- Rubber Feet (4 pieces per rack)

1. Install 4 rubber feet in the corner holes on the bottom of rack by inserting the stems through the bottom of the frame **(A)**. Pull stems until rubber feet pop into place. Excess stem may be trimmed with scissors.
2. Place grid between top and middle shelf of rack with large tabs toward center of rack **(B)**.
3. Insert large tabs onto center rib of frame **(C)**. Compress spring latch (3 small tabs) toward center of rack allowing rack to drop into outer rib **(D)**.
4. Place second grid above top shelf of frame and repeat step 3.
5. Repeat process for opposite side of rack.




Interchangeable/Replacement Grids sold separately in sets of 4

CATALOG No.	TUBE SIZE	GRID COLOR
F18745-1300	10-13mm	Orange
F18745-1600	13-16mm	Blue
F18745-2000	16-20mm	Green
F18745-2500	20-25mm	Fuchsia
F18745-3000	25-30mm	Yellow



Help others choose the right product. Write a review at www.belart.com



You May Also be Interested in:

Scienceware® Test Tube Capper

Reduce fatigue and soreness of hands and fingers caused by the repetitive task of applying test tube caps. Test Tube Capper fits comfortably in the palm of your hand. Use the tip to apply pressure and securely cap test tubes.

Catalog No. F18799-0000

661 Route 23 South, Wayne, NJ 07470

TEL: 1-800-423-5278 • FAX: 973-694-7199 • www.belart.com

Bel-Art – SP Scienceware assumes no obligation or liability for any advice furnished by it, or for results obtained with respect to these products. All such advice is given and accepted at the buyer's risk.

© 2019. All Rights Reserved.

Bel-Art® – SP Scienceware® is a division of SP Industries, Inc.



918745075 – Rev 2

