

Lab-Aire II Dryer FAQ's

Where can I mount my Lab-Aire II Dryer?

Wall-mount models have molded in tabs with key hole openings for easy hanging using standard hardware (not included) on walls or pegboard.

Does each model have pegs for every slot?

Yes, each model has a peg to go in every slot. With non-electric models you only put pegs where you need them. To maximize performance of electric models with circulating air, if you are not using a peg, you can replace a peg with a plug to block air from that slot.

What voltage are the electric Lab-Aire II Dryers available in?

Electric Lab-Aire II Dryers are available in 8-amp @120V and 4-amp at 230V. 230V Lab Aire II Dryers are CE approved.

How much weight can the pegs hold?

Small pegs (3") can hold up to 2.5 lbs.
Medium pegs (4") can hold up to 10 lbs.
Large pegs (5") can hold up to 15 lbs.

Can baskets be used on all models?

No, baskets only fit electric models which have built-in peg hooks and rests that support the baskets.

Do the Lab-Aire II Dryers have drip trays?

Yes, there is a 3" wide drip tray on each side with scalloped wells that can also be used to hold pipettes and other long glassware.

