




Lab-Aire® II Dryers Comparison Chart

Lab-Aire® II Drying Racks

Laboratory space is always in high demand making the small footprint of the Lab-Aire® II Drying Racks a practical solution for drying and storing fragile, frequently used labware.

Choose from non-electric benchtop or wallmount drying racks or accelerate drying with the warm, flowing air from electric benchtop models. All styles have removable pegs in 3 different lengths that can be positioned for optimal use. Extra accessories and a variety of different sizes/configurations let you create a drying system that works for you.

Wallmount - Non-electric	Benchtop- Non-electric	Benchtop- Electric
Space-saving, move drying off the bench	Portable, use anywhere convenience	Accelerate drying with warm, circulating air
Easily mounts on walls or standard pegboard	Offers the most choices in size	Use with or without power cycle
Single-sided; can be custom made with additional height	Single- and Double-sided models	Single- and Double-sided models: 120V or 230V



Wallmount - Non-electric
Single-sided, 2 tier



Benchtop - Non-electric
Double-sided, 1 tier



Benchtop - Electric
Single-sided, 3 tier

Increase Drying Space Three Ways

1. Go up in height – 1, 2 and 3 tier models available in all styles
2. Double-sided benchtop models double your capacity with minimal increase in footprint
3. Use multiple Lab Aire® II Dryers

Increase Flexibility

Each Lab Aire® II Dryer comes with an assortment of pegs in various lengths. Electric models also have baskets for holding small labware. Additional pegs and baskets* to customize your dryer can be purchased separately.

* Baskets fit electric models only.

Lab-Aire® II Dryers Quick Comparison Chart

Select from 12 different models with up to 114 peg positions.

- Rugged ABS plastic will never corrode and easily wipes clean
- Drip/storage trays with scalloped wells hold long thin labware
- Small (3"), Medium (4"), and Large (5") pegs can be moved and repositioned to hold a variety of labware including beakers, flasks, bottles, test tubes and cylinders
- Extra pegs and drying baskets available

Style		Wallmount	Benchtop	Electric	Dimensions (W x D x H)	Peg Slots	Peg Sizes*	Model
Wallmount Non-Electric	Single 1 Tier	X	X		37.5 x 12.7 x 24cm (14.75 x 5 x 9.4")	19	S - 6, M - 7, L - 6	F18933-0011
	Single 2 Tiers	X			37.5 x 12.7 x 41.7cm (14.75 x 5 x 16.4")	38	S - 12, M - 14, L - 12	F18933-0014
	Single 3 Tiers	X			37.5 x 12.7 x 59.4cm (14.75 x 5 x 23.4")	57	S - 18, M - 21, L - 18	F18933-0013
Benchtop Non-Electric	Single 1 Tier	X	X		37.5 x 12.7 x 24cm (14.75 x 5 x 9.4")	19	S - 6, M - 7, L - 6	F18933-0011
	Single 2 Tiers		X		37.5 x 17.8 x 39cm (14.75 x 7 x 15.4")	38	S - 12, M - 14, L - 12	F18933-0012
	Single 3 Tiers		X		37.5 x 17.8 x 57cm (14.75 x 7 x 22.4")	57	S - 18, M - 21, L - 18	F18933-0015
	Double 1 Tier		X		37.5 x 25.4 x 21.3cm (14.75 x 10 x 8.4")	38	S - 12, M - 14, L - 12	F18933-0021
	Double 2 Tiers		X		37.5 x 25.4 x 39cm (14.75 x 10 x 15.4")	76	S - 24, M - 28, L - 24	F18933-0022
	Double 3 Tiers		X		37.5 x 25.4 x 57cm (14.75 x 10 x 22.4")	114	S - 36, M - 42, L - 36	F18933-0023
Benchtop Electric	Single 2 Tiers		X	X	42.5 x 19 x 40cm (16.75 x 7.5 x 15.7")	19	S - 6, M - 7, L - 6 2 Baskets, 7 Peg Plugs	F18819-0012 (120V) F18819-2301 (230V)
	Single 3 Tiers		X	X	42.5 x 19 x 57.7cm (16.75 x 7.5 x 22.7")	38	S - 12, M - 14, L - 12 2 Baskets, 14 Peg Plugs	F18819-0013 (120V) F18819-2302 (230V)
	Double 2 Tiers		X	X	42.5 x 25.4 x 40cm (16.75 x 10 x 15.7")	38	S - 12, M - 14, L - 12 4 Baskets, 14 Peg Plugs	F18819-0022 (120V) F18819-2303 (230V)
	Double 3 Tiers		X	X	42.5 x 25.4 x 57.7cm (16.75 x 10 x 22.7")	76	S - 24, M - 28, L - 24 4 Baskets, 28 Peg Plugs	F18819-0023 (120V) F18819-2304 (230V)

© Copyright 2012 Bel-Art – SP Scienceware, Inc. All rights reserved. This information is presented in good faith. However, no warranty of any kind is made with respect to such information nor are any results guaranteed. Always read instructions for the products you are working with.

*S - Small • M - Medium • L - Large