




# Chromatographic Column Holder

Catalog Number F18202-0000

Mount the holder on two support stands by sliding the support blocks onto the support stand rods and tightening the screws.

Loops move sideways in the base and can hold columns in between loops.

When using larger diameter columns, some loops may be removed to provide more space. This should be done during set-up to prevent dislodging of hanging columns.



**BEL-ART PRODUCTS**

661 Route 23 South, Wayne, NJ 07470-6814 USA • TEL: 1-800-4BEL-ART • FAX: 973-694-7199 • [www.belart.com](http://www.belart.com)