

INSTRUCTIONS



PrepSafe™ FLOATING MICROCENTRIFUGE TUBE RACK

Catalog No. F18742-0001



INSERTING MICROCENTRIFUGE TUBES:

1. Place the rack on a flat surface.
2. Secure the rack from any movement with one hand.
3. Insert a microcentrifuge tube into the desired location.
4. Slide and push the tube completely under the immobilizing tab until the tube is completely secure.

REMOVING MICROCENTRIFUGE TUBES:

1. Place the rack on a flat surface.
2. Secure the rack from any movement with one hand.
3. Slide and push the microcentrifuge tube rack in the opposite direction until it clears the locking clips.

Bel-Art Products

661 Route 23 South, Wayne, NJ 07470-6814 USA

TEL: 1-800-4BEL-ART • FAX: 973-694-7199 • www.belart.com

918742002

Scienceware® is a registered TM of Bel-Art Products. PrepSafe™ is a TM of Bel-Art Products. Bel-Art Products assumes no obligation or liability for any advice furnished by it, or for results obtained with respect to these products. All such advice is given and accepted at the buyer's risk. © 2008 Bel-Art Products. All Rights Reserved.

12/08

INSTRUCTIONS



PrepSafe™ FLOATING MICROCENTRIFUGE TUBE RACK

Catalog No. F18742-0001



INSERTING MICROCENTRIFUGE TUBES:

1. Place the rack on a flat surface.
2. Secure the rack from any movement with one hand.
3. Insert a microcentrifuge tube into the desired location.
4. Slide and push the tube completely under the immobilizing tab until the tube is completely secure.

REMOVING MICROCENTRIFUGE TUBES:

1. Place the rack on a flat surface.
2. Secure the rack from any movement with one hand.
3. Slide and push the microcentrifuge tube rack in the opposite direction until it clears the locking clips.

Bel-Art Products

661 Route 23 South, Wayne, NJ 07470-6814 USA

TEL: 1-800-4BEL-ART • FAX: 973-694-7199 • www.belart.com

918742002

Scienceware® is a registered TM of Bel-Art Products. PrepSafe™ is a TM of Bel-Art Products. Bel-Art Products assumes no obligation or liability for any advice furnished by it, or for results obtained with respect to these products. All such advice is given and accepted at the buyer's risk. © 2008 Bel-Art Products. All Rights Reserved.

12/08