

**WARNING
PLEASE NOTE!**



Safety Acid/Solvent Bottle Carriers

**Catalog No. 16957-0000, 16958-0000, 16959-0000,
16960-0000**

Scienceware® Safety Acid/Solvent Bottle Carriers are made to hold either a standard U.S. 5 pint reagent bottle with maximum diameter of 5" (127 mm) or a standard U.S 1 gallon reagent bottle with maximum diameter of 6" (159 mm).

Certain bottles of non-U.S.A. manufacturers will not fit the standard Scienceware® Safety Acid/Solvent Bottle Carriers.

BEL-ART PRODUCTS

WAYNE, NJ 07470 USA

TEL: 1-800-4BEL-ART • FAX: 973 694-7199 • www.belart.com

962150028

11/02

Scienceware®, is a reg. TM of Bel-Art Products
Bel-Art Products assumes no obligation or liability for any advice furnished by it, or for results obtained with respect to these products. All such advice is given and accepted at the buyer's risk.

**OTHER FINE SCIENCEWARE® PRODUCTS
FOR YOUR LABORATORY**

Eye Wash Station



Double eyewash safety station pictured, single station also available.

Bright, Easy to Locate

These easily recognized, ready-to-use Eye Wash Stations come complete with eye wash bottles, dust covers, mirror, and instructions for use and service. A bracket holds the bottles securely in place, yet bottles are easily released when pulled out. Includes station and empty eye wash bottles. We recommend filling with sterile eyewash solution catalog number 24879-0032.

Catalog No.

24865-0000
24866-0000
24868-0000

Bottle Size

Station with one 500ml (16oz) empty eye wash bottle
Station with one 1000ml (32oz) empty eye wash bottle
Station with two 1000ml (32oz) empty eye wash bottle

Recommended Refill

24879-0032

(32oz) sterile eyewash solution refill



BEL-ART PRODUCTS

661 Route 23 South, Wayne, NJ 07470 USA

Tel: 1-800-4BELART • Fax: 973-694-7199 • www.belart.com