



Safety Step Stool

Catalog No. F24900-0000



The Safety Step Stool features a wrap-around rubber bottom skirt that gives a non-skid contact with the floor. The stool rolls on three spring-loaded swiveling casters. When weight is placed on the stool, the casters recede against compression springs and the stool will rest directly on the floor. When weight is removed, the springs push the casters out for easy rolling on the floor. The stool rolls easily on any smooth hard surface and will also roll well on low-pile carpets.

The Safety Step Stool conforms to EN 14183:2003 testing standard and is approved for a weight bearing capacity of 330lb (150kg). The stool is molded of rugged ABS plastic and can be cleaned with water and detergent or soap. Do not use petroleum-based solvents.

BEL-ART PRODUCTS

Wayne, NJ 07470 USA • TEL: 1-800-4BEL-ART • FAX: 973-694-7199 • www.belart.com

999320738

© 2014 Bel-Art Products. All rights reserved. Scienceware®, is a reg. TM of Bel-Art Products

Bel-Art Products assumes no obligation or liability for any advice furnished by it, or for results obtained with respect to these products. All such advice is given and accepted at the buyer's risk.

01/14

