

INSTRUCTIONS



Eye Wash Bottle Wall Rack

Catalog No. 24859-0000



The epoxy coated steel rack holds a single 32 ounce (1000cc) Eye Wash Bottle, F24851-0000. It is designed to be hung on a wall by means of a single screw or nail. It is important that the screw be anchored properly, into a wooden stud or into a plastic or lead anchor if no stud is available. If a nail is used it must go into a solid wooden member. Other anchoring systems may be used for wallboard or plaster walls.

BEL-ART PRODUCTS

661 Route 23 South, Wayne, NJ 07470-6814 USA
TEL: 1-800-4BEL-ART • FAX: 973-694-7199 • www.belart.com

999320809

07/02

Scienceware®, is a reg. TM of Bel-Art Products
Bel-Art Products assumes no obligation or liability for any advice furnished by it, or for results obtained with respect to these products.
All such advice is given and accepted at the buyer's risk.

OTHER FINE SCIENCEWARE® PRODUCTS FOR YOUR LABORATORY

Eye Wash Station

Catalog No. 24865-0000, 24866-0000



Bright, Easy to Locate

This easily recognized, ready-to-use Eye Wash Station comes complete with eye wash bottle, dust cover, mirror, and instructions for use and service. A bracket holds the bottle securely in place, yet the bottle is easily released when pulled out. Backboard is 44 x 38cm (17½" x 15") with 4 screw holes for wall mounting. Includes station and bottle only, eye wash not included. Bottles supplied empty.

Catalog No.
F24865-0000
F24866-0000

Bottle Size
Station with 500ml (16oz) empty eye wash bottle
Station with 1000ml (32oz) empty eye wash bottle



BEL-ART PRODUCTS

661 Route 23 South, Wayne, NJ 07470-6814 USA
TEL: 1-800-4BEL-ART • FAX: 973-694-7199 • www.belart.com