

Press Release
9/16/2015

For Immediate Release
Contact:
Teri L. Rocco
trocco@belart.com – 1-800-423-5278
Bel-Art – SP Scienceware

New Freezer Boxes Save Big on Storage Space and Shipping Costs

Wayne, NJ – September 16, 2015 – Storing 100 empty cardboard freezer boxes is a challenge even for a large size laboratory, not to mention the additional cost of shipping a large package essentially filled with air. New Pop Up 2” Freezer Boxes from Bel-Art – SP Scienceware ship flat and store flat when not in use, reducing your needed storage space and saving you money on shipping costs.

“These boxes assemble in seconds and are the same size as industry standard 2” cardboard freezer boxes for easy integration into the laboratory,” according Mike Grzelak, Product Manager for Bel-Art – SP Scienceware. “However, unlike standard cardboard freezer boxes, these boxes are more durable and moisture-proof, so they will not warp, bend, or even stick to icy freezer racks. This helps to save even more money as the need to replace freezer boxes will be less often.”

Made of translucent polypropylene, the Pop Up 2” Freezer Boxes are resistant to mold and will withstand temperatures down to -86°C. Use the box with or without the optional 81-place divider that holds 2ml cryo-tubes and has a numerical grid for easy sample identification.

Suggested list price for the Scienceware® Pop Up 2” Freezer Box with 81-place divider, catalog no. F18829-0012, is



661 Route 23 South, Wayne, NJ 07470
Tel: 973-694-0500 • Fax: 973-694-7199 • email: info@belart.com • www.belart.com



\$62.50 for a pack of 12. Boxes and dividers are also sold separately. The Pop Up Box, catalog no. F18829-0000, is great for storing more than just vials and is \$30.00 for a pack of 12. The 81-Place Divider, catalog no. F18829-0001, fits any 5¼" x 5¼" freezer box and is \$35.00 for a pack of 12.

To see how easily the Scienceware® Pop Up 2" Freezer Boxes assemble, visit the [Scienceware YouTube channel](http://www.belart.com/popup) (www.belart.com/popup).

For more information on the Scienceware® Pop Up 2" Freezer Boxes, or the complete line of laboratory offerings from Bel-Art – SP Scienceware, call 1-800-423-5278, email info@belart.com or visit www.belart.com. To download catalogs from any of the SP Scienceware brands (Bel-Art, H-B Instrument or Wilmad-LabGlass), visit www.belart.com/catalogs.

Research, industrial and healthcare professionals rely on SP Scienceware for their sample handling needs. Sold worldwide, Bel-Art and H-B Instrument brand laboratory supplies include over 5,000 items from safety wash bottles to fume hoods, thermometers, racks, desiccators, hydrometers, ice buckets, scoops and spoons. Custom calibration services, plastic molding, specialty coatings and wire manufacturing are also available. Wilmad-LabGlass brand laboratory glassware complements Bel-Art and H-B Instrument with over 7000 items including NMR and EPR consumables and specialty glassware, as well as custom fabrication for individual glassware, glass repair services and OEM glass parts.

SP Scienceware is a division of Warminster, PA based SP Industries which also include SP Scientific laboratory, thermal management, and pharmaceutical production equipment and SP Ableware aids for daily living.

#

661 Route 23 South, Wayne, NJ 07470
Tel: 973-694-0500 • Fax: 973-694-7199 • email: info@belart.com • www.belart.com

An SP Industries Brand