

Press Release
August 2, 2016

For Immediate Release
Contact:
Teri L. Rocco
trocco@belart.com – 1-800-423-5278
Bel-Art – SP Scienceware

Read Temperature and Humidity Data on Your Mobile Device with the new DURAC® Bluetooth Thermometer Hygrometer

Wayne, NJ – August 2, 2016 – Remote monitoring and logging of critical environmental data in rooms or chambers is now possible with the H-B Instrument [DURAC® Bluetooth Thermometer - Hygrometers](#). These new units record up to 30 days of temperature and humidity, or temperature data only and transmit the information wirelessly, eliminating the need for a physical connection.

“Readings can be broadcast to any apple or android smartphone or device with Bluetooth 4.0 or better at distances up to approximately 40 meters utilizing the free THERMSmart app,” explained Staci Hiergesell, Product Manager for SP Scienceware.



Hourly data is recorded for up to seven days, then maximum and minimum temperatures (model B60911-0100) or temperature and humidity data (model B61506-5100) are recorded for the next 23 days. An alarm will sound if the temperature falls outside of the set parameters.

“Data is viewable as a chart or graph, or on a calendar, and data from up to six data loggers can be plotted on one graph for direct comparison,” added Hiergesell. “Data can also be downloaded as a CSV file if needed.”

661 Route 23 South, Wayne, NJ 07470
Tel: 973-694-0500 • Fax: 973-694-7199 • email: info@belart.com • www.belart.com



Both models include a Certificate of Compliance indicating accuracy for 4 temperature points or 4 temperature and 3 humidity points. They are individually serialized, CE marked and RoHS compliant.

The Bluetooth Thermometer Hygrometer, catalog no. [B61506-5100](#), has a temperature range of -10/50°C (14/122°F) and a humidity range of 10/99%, and a suggested list price of \$110.96. The Bluetooth Thermometer, catalog no. [B60911-0100](#), has a range of -10/50°C (14/122°F) and a suggested list price of \$101.96.

For more information on H-B Instrument [DURAC® Bluetooth Thermometer Hygrometers](#), or the complete line of laboratory offerings from Bel-Art - SP Scienceware, call 1-800-423-5278, email info@belart.com or visit www.belart.com. To download catalogs for any of the SP Scienceware brands (Bel-Art, H-B Instrument, or Wilmad-LabGlass), visit www.belart.com/catalogs.

Research, industrial and healthcare professionals rely on SP Scienceware for their sample handling needs. Sold worldwide, Bel-Art and H-B Instrument brand laboratory supplies include over 5,000 items from safety wash bottles to fume hoods, thermometers, racks, desiccators, hydrometers, ice buckets, scoops and spoons. Custom calibration services, plastic molding, specialty coatings and wire manufacturing are also available. Wilmad-LabGlass brand laboratory glassware complements Bel-Art and H-B Instrument with over 7000 items including NMR and EPR consumables and specialty glassware, as well as custom fabrication for individual glassware, glass repair services and OEM glass parts.

SP Scienceware is a division of Warminster, PA based SP Industries which also include SP Scientific laboratory, thermal management, and pharmaceutical production equipment and SP Ableware aids for daily living.

#

The Bluetooth® word mark and logos are registered trademarks owned by the Bluetooth SIG, Inc. and any use of such marks is under license. Other trademarks and trade names indicated are property of Bel-Art – SP Scienceware.

661 Route 23 South, Wayne, NJ 07470
Tel: 973-694-0500 • Fax: 973-694-7199 • email: info@belart.com • www.belart.com

An SP Industries Brand