

INSTRUCTIONS



Fluorescent Ruler

Catalog Number 13355-0001



The Fluorescent ruler can be in observing and photographing Ethidium Bromide stained nucleic acid agarose and acrylamide gels. The scale will fluoresce visibly upon illumination with ultraviolet light. This serves as a measure of nucleic acid mobility, and can be used to quantify the sizes of the nucleic acid fragments. This measurement will also allow the comparison of gels run at different times under the same conditions.

BEL-ART PRODUCTS

661 Route 23 South, Wayne, NJ 07470-6814 USA
TEL: 1-800-4BEL-ART • FAX: 973 694-7199 • www.belart.com

913355003

02/12

Scienceware® is a reg. TM of Bel-Art Products
Bel-Art Products assumes no obligation or liability for any advice furnished by it, or for results obtained with respect to these products. All such advice is given and accepted at the buyer's risk.