




Rotary Pipette Stand Catalog No. 18957-0000



Using the photograph as a guide, press the disc with concentric grooves onto one end of the short larger diameter tube, then press the disc with multiple holes onto the other end of the short, large diameter tube. Make sure these three parts are firmly pressed together

Slide this disc and tube assembly over the longer small diameter tube.

Press the open end of the longer tube with disc assembly into the base center.

All parts of this rack are made of polypropylene and autoclavable

661 Route 23 South, Wayne, NJ 07470



TEL: 1-800-423-5278 • FAX: 973-694-7199 • www.belart.com

918957001 - 05/15

Bel-Art - SP Scienceware assumes no obligation or liability for any advice furnished by it, or for results obtained with respect to these products.

All such advice is given and accepted at the buyer's risk.

© 2015. All Rights Reserved.

Bel-Art® - SP Scienceware® is a division of SP Industries, Inc.

