


INSTRUCTIONS



Velveteen Squares

Catalog No. H37848-0001 (Off White) & H37848-0002 (Black)



These 100% cotton, sterile Velveteen Replica Plating Squares are especially suited for use in conjunction with the Scienceware® Replica Plating Tool #F37848-0000.

The Velveteen Squares are supplied sterilized and ready for use. These may be re-used a number of times after careful washing and sterilization. To maintain sterility do not open the remaining bags until they are ready to be used. Seal tightly when not in use to reduce airborne contamination.



661 Route 23 South, Wayne, NJ 07470-6814 USA
TEL: 1-800-4BEL-ART • FAX: 973-694-7199 • www.belart.com

937848016
Part of SP Industries
Bel-Art assumes no obligation or liability for any advice furnished by it, or for results obtained with respect to these products. All such advice is given and accepted at the buyer's risk.

04/15

OTHER FINE SCIENCEWARE® PRODUCTS FOR YOUR LABORATORY

Replicate! Replicate! Replicate!

The Bel-Blotter™ 96 Wells Replicating Tool

Catalog No. F37876-0002



- **Efficient, Cost Effective, Precise**
- **Capillary Action Insures That Each Well Will Pick Up 10µl Per Tip**
- **96 "Open Ended" Pipette Tips Mate Perfectly With Any 96 Well Plate**
- **Molded of Autoclaveable Polypropylene**



661 Route 23 South, Wayne, NJ 07470-6814 USA
TEL: 1-800-4BEL-ART • FAX: 973-694-7199 • www.belart.com