

INSTRUCTIONS



Flexible Arm With Accessories

Catalog Number F18315-0000 & F18317-0000



The Flexible Arm (available separately, catalog number F18315-0000) has a ½" stainless steel pin at one end for mounting into the weighted base. Insert the end with stainless steel pin into any of the 4 wells in the weighted base and secure by tightening the thumbscrew. The base will hold a maximum of four of the Flexible Arms.

The accessory Electrode Holder (catalog number F18317-1000) snaps onto the end of the flexible arm and the electrode, reference electrode, and ATC probe are placed in appropriate positions.

BEL-ART PRODUCTS

661 Route 23 South, Wayne, NJ 07470-6814 USA

TEL: 1-800-4BEL-ART • FAX: 973 694-7199 • www.belart.com

999320808

Scienceware® is a reg. TM of Bel-Art Products

Bel-Art Products assumes no obligation or liability for any advice furnished by it, or for results obtained with respect to these products. All such advice is given and accepted at the buyer's risk.

03/12

OTHER FINE SCIENCEWARE® PRODUCTS FOR YOUR LABORATORY

Cleanware™ Electrode Washer

Catalog Number: 17188-0000

Save time and reduce scratching.

The Electrode Washer provides a safe, convenient way to clean delicate pH electrodes in place without the risk of scratching or damage. The unit consists of a polyethylene spray chamber to surround electrodes, a 500ml polyethylene reservoir that holds distilled water and a length of tubing to drain the chamber. The electrode hangs into the spray chamber and with a squeeze of the reservoir, a spray of water rises the electrode. The base of the system is 14cm (5 1/2") diameter and the washer is 21cm (8 1/4") high.



Electrode Rack

Catalog Number: 18509-0000

This handy polypropylene, autoclavable rack holds 8 test tubes for storing and conditioning pH and ion electrodes. Holds 2 rows of 4 tubes up to 20mm in diameter.



BEL-ART PRODUCTS

661 Route 23 South, Wayne, NJ 07470-6814 USA

TEL: 1-800-4BEL-ART • FAX: 973 694-7199 • www.belart.com