




Klik-Klip™ Wire and Tubing Clip
 Catalog No. F18234-0000



1. Determine the exact location to place the wire clip.
 2. Surface should be clean, dry, and wax-free
 3. Peel off backing paper and press the wire clip onto the surface. Install with open end up. Press firmly for 30 seconds.
 4. Allow 2 hours for the adhesive to set.
- NOTE: Once in place, the clip cannot be moved without damage to the adhesive backing.**
5. Open the clip and put the wire or tube in place, then close the clip.

Help others choose the right product. Write a review at www.belart.com




661 Route 23 South, Wayne, NJ 07470
 TEL: 973-694-0500 • FAX: 973-694-7199 • www.belart.com

918234001 – Rev 2

Klik-Klip™ is a trademark of Bel-Art – SP Scienceware
 Bel-Art – SP Scienceware assumes no obligation or liability for any advice furnished by it, or for results obtained with respect to these products.
 All such advice is given and accepted at the buyer's risk.
 © 2018. All Rights Reserved.
 Bel-Art® – SP Scienceware® is a division of SP Industries, Inc.



Klik-Klip™ Wire and Tubing Clip
 Catalog No. F18234-0000



1. Determine the exact location to place the wire clip.
 2. Surface should be clean, dry, and wax-free
 3. Peel off backing paper and press the wire clip onto the surface. Install with open end up. Press firmly for 30 seconds.
 4. Allow 2 hours for the adhesive to set.
- NOTE: Once in place, the clip cannot be moved without damage to the adhesive backing.**
5. Open the clip and put the wire or tube in place, then close the clip.

Help others choose the right product. Write a review at www.belart.com



661 Route 23 South, Wayne, NJ 07470
 TEL: 973-694-0500 • FAX: 973-694-7199 • www.belart.com

918234001 – Rev 2

Klik-Klip™ is a trademark of Bel-Art – SP Scienceware
 Bel-Art – SP Scienceware assumes no obligation or liability for any advice furnished by it, or for results obtained with respect to these products.
 All such advice is given and accepted at the buyer's risk.
 © 2018. All Rights Reserved.
 Bel-Art® – SP Scienceware® is a division of SP Industries, Inc.

