


INSTRUCTIONS



Lab-Aire® II BENCHTOP LABWARE DRYING RACK

Catalog No. 18933-0012 • 18933-0015 • 18933-0021 • 18933-0022 • 18933-0023



PEG ASSEMBLIES:

On each panel, pegs should be arranged with the small pegs on the top rows and large pegs on the bottom rows. To ensure that pegs are not impacted too tightly for removal or rearrangement, do not force the pegs into the peg holes.

REPLACEMENT PARTS

<u>DESCRIPTION</u>	<u>PART NO.</u>
Bag of 19 Small Pegs & 6 Plugs	18931-0001
Bag of 19 Medium Pegs & 6 Plugs	18931-0002
Bag of 19 Large Pegs & 6 Plugs	18931-0003



661 Route 23 South, Wayne, NJ 07470
TEL: 1-800-4BEL-ART • FAX: 973-694-7199 • www.belart.com

918927115

Bel-Art – SP Scienceware assumes no obligation or liability for any advice furnished by it, or for results obtained with respect to these products. All such advice is given and accepted at the buyer's risk.
© 2015. All Rights Reserved.

Bel-Art® – SP Scienceware® is a division of SP Industries, Inc.

05/15