

I N S T R U C T I O N S




60mm Petri Dish Rack

Catalog No.
F18991-0060

ASSEMBLY:

1. Refer to drawing for orientation of parts.
2. Assemble the center divider to the base securing with one nut on each of the two threaded posts.
3. Assemble each of the rods to the base, secure each rod with one nut on the threaded post.



BEL-ART PRODUCTS

661 Route 23 South, Wayne, NJ 07470-6814 USA
TEL: 1-800-4BEL-ART • FAX: 973-694-7199 • www.belart.com

OTHER FINE SCIENCEWARE® PRODUCTS FOR YOUR LABORATORY

Write-On™ Tape Dispensers make labeling easy and fast


Producing legible labels is so much easier on a smooth, flat platform than on curved surfaces such as bottles or vials. The new Scienceware® Write-On™ Tape Dispenser provides this convenience by always having labeling tape positioned on the large writing platform, ready to make into a label. Select either the benchtop model with weighted base or economy model that can be wall mounted to save space. All models are easy to use by both left and right handed users.

Each dispenser can hold 5 or more tape rolls with widths from 12 - 25mm (½ - 1"). Ideal for use with the full line of Scienceware, Write-On™ Tapes, F13463 series, available in 6 colors and three widths, 12, 19, and 25mm (½", ¾" and 1"). This advanced tape adheres to any clean surface and resists moisture and most solvents. It withstands long exposure to temperatures as low as -73°C (-68°F) and up to 80°C (176°F) yet peels off without leaving a sticky residue.

Contact your preferred distributor to order these and other Scienceware® products.

CATALOG NO.	DESCRIPTION	DIMENSION
F13472-0000	Benchtop Model	302 x 143 x 121 mm H (12 x 5 ⁵ / ₈ x 4 ³ / ₁₆ " H)
F13473-0000	Economy Model	302 x 143 x 106 mm H (12 x 5 ⁵ / ₈ x 4 ³ / ₁₆ " H)

Economy Model (shown)
Can be wall-mounted



Benchtop Model (not shown)
Includes weighted base



BEL-ART PRODUCTS

661 Route 23 South, Wayne, NJ 07470-6814 USA
TEL: 1-800-4BEL-ART • FAX: 973-694-7199 • www.belart.com

918991066

Scienceware®, is a reg. TM of Bel-Art Products
Bel-Art Products assumes no obligation or liability for any advice furnished by it, or for results obtained with respect to these products. All such advice is given and accepted at the buyer's risk.

03/09