


INSTRUCTIONS

Buchner Funnel Trowel

Catalog No. F36820-0118



The Bel-Art – SP Scienceware® Buchner Funnel Trowel is designed to smooth the cake (i.e. the precipitate) in 18, 24 and 36" Buchner Funnels, allowing for more efficient filtration.

Directions for Use:

Begin filtration in the Buchner funnel. Holding the Bel-Art – SP Scienceware® Buchner Funnel Trowel at a slight angle in either left or right hand, sweep the trowel over the surface of the cake and smooth over any noticeable cracks or fissures. Continue smoothing the cake as needed until desired filtration level is reached.

Care and Cleaning:

- Do NOT autoclave product.
- Compatible with industrial dishwasher water temperature and cleaners
- Recommended detergent: Aquet®, Catalog no. F17094-0030, 1 gallon

661 Route 23 South, Wayne, NJ 07470

TEL: 973-694-0500 • FAX: 973-694-7199 • www.belart.com

Bel-Art – SP Scienceware assumes no obligation or liability for any advice furnished by it, or for results obtained with respect to these products.

All such advice is given and accepted at the buyer's risk.

© 2015. All Rights Reserved.

Bel-Art® – SP Scienceware® is a division of SP Industries, Inc.




936820119 – 10/15



INSTRUCTIONS

Buchner Funnel Trowel

Catalog No. F36820-0118



The Bel-Art – SP Scienceware® Buchner Funnel Trowel is designed to smooth the cake (i.e. the precipitate) in 18, 24 and 36" Buchner Funnels, allowing for more efficient filtration.

Directions for Use:

Begin filtration in the Buchner funnel. Holding the Bel-Art – SP Scienceware® Buchner Funnel Trowel at a slight angle in either left or right hand, sweep the trowel over the surface of the cake and smooth over any noticeable cracks or fissures. Continue smoothing the cake as needed until desired filtration level is reached.

Care and Cleaning:

- Do NOT autoclave product.
- Compatible with industrial dishwasher water temperature and cleaners
- Recommended detergent: Aquet®, Catalog no. F17094-0030, 1 gallon

661 Route 23 South, Wayne, NJ 07470

TEL: 973-694-0500 • FAX: 973-694-7199 • www.belart.com

Bel-Art – SP Scienceware assumes no obligation or liability for any advice furnished by it, or for results obtained with respect to these products.

All such advice is given and accepted at the buyer's risk.

© 2015. All Rights Reserved.

Bel-Art® – SP Scienceware® is a division of SP Industries, Inc.



936820119 – 10/15

