


Mini LED Light Box

Catalog Number F37864-3000

- Remove battery cover
- Insert 6 AAA batteries (not included); can also be powered by 120V AC adaptor (not included)
- Use switch on lower left side to turn box on
- Gently place grid to cover entirety of viewing area (if applicable)



You May Also be Interested in:



Extra Counting Grids (3)

F37864-3200



Colony Counter System

F37863-3000

complete system with magnifier, colony counter, and light box



Mini Magnifier

F37865-0000

for further enhanced viewing



Colony Counter

F37862-0002

for keeping track of colonies digitally

Your opinion counts! Write a product review at www.sp-belart.com



Bel-Art Products, a part of SP

661 Route 23 South, Wayne, NJ 07470 USA

Tel: 973-694-0500 | cservice.belart@spindustries.com

sp-belart.com

Individual use varies. Bel-Art Products/SP assumes no obligation or liability for any advice furnished by it or for results obtained with respect to these products. All advice is given and accepted at the user's risk. For warranty information visit belart.com/corporate.

