


# I N S T R U C T I O N S



## Spinbar® Magnetic Sink Strainer

Catalog No. F37787-0000 (blue) • F37787-0001 (green)



### INSTRUCTIONS:

- Place strainer in sink drain in an upright position so that Scienceware® logo is visible.
- Flat edge allows strainer to be placed close to sink walls.
- Four strong magnets in each corner attract loose Spinbar® magnetic stirring bars and prevent them from falling into the drain.
- Periodically inspect and clean strainer as needed.



## BEL-ART PRODUCTS

PEQUANNOCK, NJ 07440-1992 USA

TEL: 1 800 4BEL-ART • FAX: 973 694-7199 • [www.belart.com](http://www.belart.com)

937902113

01/08

Scienceware® is a registered TM of Bel-Art Products • Spinbar® is a registered TM of Bel-Art Products  
Bel-Art Products assumes no obligation or liability for any advice furnished by it, or for results obtained with respect to these products. All such advice is given and accepted at the buyer's risk. © 2007 Bel-Art Products. All Rights Reserved.