

Application Focus

Metal detectable plastic scoops help prevent recalls and operational delays

Challenge:

Sampling tools and plastic fragments are left behind in bags and containers

Application:

Quality-Control

Sampling products can be forgotten and left inside containers. As a result, finished products packages can contain sampling scoops and/or plastic fragments, leading to expensive recalls and operational delays.



Solution:

Sterileware® Sense-able Blue Metal Detectable Scoops H36902-0006, H36904-0006, H36906-0006

The Sterileware Sense-able scoops are detectable by standard X-ray and metal detectors. Run any raw good or finished product packages through X-ray/Metal detectors and quickly find any sampling scoops or fragments that were accidentally left behind.

Polystyrene plastic impregnated with stainless steel powder from top to bottom permits scoops and fragments as small as 1.6mm to be detected. Plants can now QC check with optimal time and efficiency.

For more information visit:

<https://www.belart.com/sterileware-sense-able-scoops-metal-detectable-sterile-sampling-tools.html>

